

QUALITY
MADE IN GERMANY



New: Universal Smartphone Adapter mag with two different fixing buttons

Dear HEINE partner,

effective immediately, the [Universal Smartphone Adapter mag](#) for dermatoscopes HEINE DELTA 30 PRO | DELTA 30 and DELTAone is delivered with a set of two fixing buttons of different heights. Users can use almost any smartphone or tablet with the adapter, regardless of the camera position.

The low fixing button is suitable for most models. The higher fixing button is better suitable for newer models with higher camera rims, e.g. the iPhone 16 PRO.



K-000.34.273

The adapter is delivered with two fixing buttons now.

K-000.34.274

The fixing buttons are of different depth: lower and higher button (L). The set can be ordered separately if needed.

Need a new fixing button? Simply order the set separately.

The set that contains the two different fixing buttons can be ordered separately from the adapter. The catalogue number remains the same as before: K-000.34.274.

Check out our instructional [video how to fix a button to the smartphone and how to adjust the adapter](#).

Please remember to update your marketing materials accordingly.
We are happy to answer any other questions you may have.

Best wishes from Gilching
Team HEINE Optotechnik



HEINE Optotechnik
Dornierstraße 6, 82205 Gilching, Germany